


# EZFeed

Peg Kit



An ISO 9001-2008 Company



The PEG KIT is placed into the stomach percutaneously for nutrition and medications in patients who cannot swallow.



# EZFeed

## Peg Kit

The peg tube

helps in the transport of food and other liquids

The feeding adapter

helps for the inlet of medicine and other fluids into the peg tube

The bolster

is used to affix the tube firmly

The blue pull wire, scalpel and introducer needle

helps in the procedural set up of peg tube

## Product Specifications

<b>Ref No.</b>	702100
<b>Item</b>	PEG SET 24F
<b>Contents</b>	1 no each of peg tube, feeding adapter, bolster, scalpel, clamp, introducer needle and blue pull wire.
<b>Ref No.</b>	701100
<b>Item</b>	PEG SET 20F
<b>Contents</b>	1 no each of peg tube, feeding adapter, bolster, scalpel, clamp, introducer needle and blue pull wire.

Regd. Office - 344/5, N.S.C. Bose Road, Kolkata - 700 047, India. Tel - (033)2421 5614

Website : [www.imi.co.in](http://www.imi.co.in)

European Authorised Representative : Medical Technology Promedt Consulting GmbH

Altenhofstrasse 80 D-66386 St. Ingbert, Germany

Tel : +49 (0) 6894 581020 Fax: +49 (0) 6894 581021



013



0434



A SUBSIDIARY OF VITAE CORE COMPANY



An ISO 9001-2008 Company